

Attention

1. Please take off shoes before taking measurement to avoid scraping or scratching the glass platform.
2. The glass platform may be pretty slippery when wet, thus please make sure both the glass and your feet are dry before using.
3. Please use soft tissue with alcohol or glass cleaner to clean the surface if it is dirty. No soap or other chemicals. Keep it away from water, heat and extreme coldness.
4. The scale is high precision measuring device. Never jump or stomp on the scale or disassemble it and please handle it carefully when moving it.
5. The scale is for family use only and not suitable for professional use.
6. Please check the battery if the scale malfunctions. Change a new piece if necessary. For other problems, please contact the local dealer or our company directly by phone or E-mails.



TONE

Talking Bathroom Scale

USER'S MANUAL



Thank you for choosing this electronic personal scale. To ensure proper operation and to maintain the lifespan of the scale, please read this manual carefully before use.

Key features:

1. High precision strain gauge sensors system.
2. Tempered-glass platform.
3. Capacity: 150kg, 180kg
4. Division: 0.1 kg
5. LCD display
6. Auto on.
7. Auto zero & Auto off.
8. Low battery and over load indication.
9. Free switch within three different units: kg/lb/st.
10. Battery function: 3*AAA batteries (included).
11. True voice data reading function

Operation

1. Find the lithium cell in the packet and battery holder at the back of the scale, open the cover, install the battery into the holder and cover it back (make sure the positive pole of the battery upwards).
2. Place the scale on hard and flat surface. Clean floor is best.
3. Step on the scale directly and the weight will display.



4. Stay still on the glass platform. During weighting, weight data will be read by 3 times, then the unit will switch off automatically.

Caution

1. To ensure accuracy, we suggest you take measurement at a fixed time every day (e.g. 7:00pm-9:00pm).
2. Please change a new battery when "LO" appears on the LCD which indicates that the battery is out of power.



3. "Err" showed on LCD indicates overload .

