TONE

INSTRUCTIONS

Electronic Talking Scale

- 3. Please stand steady while on scale to obtain the proper weight. Scale will display weight once stable.
- 4. After the weight is locked, the electronic scale will broadcast (speak) the weight according to display weight. The voice can be changed to English or Spanish. If you do not wish to hear the weight, set the voice selection to OFF.



Automatic shutdown

After weighing, the scale will automatically turn off. If scale is turned on but remains at zero, scale will turn off in approximately six seconds as an energy-saving design.

Overload indicator

If weight above the full range of the scale is reached, the scale will show the symbol below:



Basic characteristics

Product is a high precision talking electronic scale

- Includes high accuracy sensors
- Weighs 400lbs / 180kg
- Weighing units in lbs / kgs
- LCD display
- Automatic zero tracking
- Holds weight on display once stable
- Automatic shutdown
- Overload indicator

Power

This product uses low power consumption in an environmentally friendly design. Scale comes with battery. When the scale shows low power, please replace the battery.

Method of use

- 1. Weigh on a level surface.
- 2. Scale will check weight twice, then lock the weight on the display



Warnings:

- 1. The scale should be kept dry, avoiding water to prevent slipping while weighing.
- 2. Use this scale under normal weighing purposes only. Do not hit or drop scale to avoid damage or breakage.
- 3. Do not use chemical solvents on the scale. When necessary, wipe with a damp cloth.
- 4. Do not store in a damp environment
- 5. Please store safely when not needed for long term. Remove battery while scale is being stored.
- 6. If scale does not turn on immediately, apply extra pressure. If scale still will not turn on, check battery connection. If necessary change battery.
- 7. This scale is for consumer use only and should not be used in a commercial environment.